

No-needle Meso Series

No-needle Meso Therapy System

Applications

- Nutrition guide in
- Skin whitening,
- Skin tightening, lifting
- Face contouring and reshaping.
- Removing wrinkles, fine lines. Shrinking large pores.
- Improving elasticity and tone of skin.



How does it work

No-needle eMeso therapy, adopting special electric current pulse, stimulate and create Electrophoresis reactions on skin cells, which will help open protein and nutrition acceptance channels. Thus, cells activity and permeability of skin can be improved.

Nutrition can be conducted in and improve skin conditions

Main Features

- Special no-needle meso treatment handpiece with magnetic ball.
- Vacuum massage and sucking lifting, faster result getting.
- 6 exchangeable vacuum treatment heads for various positions.
- Blue touch LCD with user friendly icon menu.



Specifications

Power	AC 220~240 V, 50 Hz	Working Humidity	≤80%
Vacuum	Max. 65Kpa, adjustable	Working Noise	≤65 db
Micro-current	99 levels adjustable	Working Temperature	5 ~ 40℃
Dimensions / Weight	350 x 440 x 250 mm / 8Kg		