

Oxygen Jet Peel



THEORY


A revolutionary patented peeling system with cutting-edge technology based on modern jet aviation principle. This treatment uses only oxygen and water or saline solution micro-droplets (accelerated to supersonic speeds) to achieve the desired effect. It leaves no burning sensation or redness. It leaves the skin hydrated smooth and polished.

FUNCTION

Skin care and skin testing system.

REVOLUTIONARY APPROACH FOR

1. Sun damage and age spots
2. wrinkle reduction
3. Stretch mark reduction
4. Acne and acne scar reduction
5. Calming rosacea
6. Evens pigment
7. Unblock pores
8. Skin polishing
9. Supersonic saline skin exfoliation
10. Deeply infuse pure oxygen
11. Safe on all skin types



TECHNICAL CHARACTERISTICS

Screen	LCD
Oxygen Source	Medical oxygen class(local supply)
Related working pressure range	0.3 0.5Mpa
Consumption of saline water	1.5 to 2.5 ml/minute
Consumption of medical oxygen	25 35 L/minute
Injecting velocity	200m/s
Skin testing pixel	1. 3 mega pixels
Skin testing probe	Derma and epidermal layer interconvert
Size	
G.W	
Voltage	220V/110V, 60Hz/50Hz.
Certification	CE
OEM or ODM	Yes